

La Bella Strings

Guide: Finding the Right Scale Set for Your Bass

Knowing your bass' "scale length" is not enough information to purchase the correct set.


Figuring out your bass' required "**string wound length**" is key BEFORE you purchase a set from ANY string manufacturer.

How To Figure out the "String Wound Length" your Bass Requires:

On your bass, measure the distance from the ball-end to $\frac{1}{4}$ " past the nut (you may have to remove the string if the string goes through the body).

Then use our key to figure out the correct set for your bass:

La Bella Short Scale Set:	32" $\frac{1}{4}$ " Wound Length
La Bella Medium Scale Set:	34" $\frac{1}{4}$ " Wound Length
La Bella Standard/Long Scale Set:	37 $\frac{1}{4}$ " Wound Length
La Bella Extra Long Scale:	38" Wound Length



Note of the Following:

- Make sure the string wound length clears the nut.
- For flat wounds, make sure the wound portion of the string does NOT wrap around the tuning post (or else they will break).
- Do not use standard flat wound bass sets on basses that require stringing through the body of the bass (or else they will break). We have Through-Body "-TB" models for those basses e.g. 760FL-TB, 760FS-B-TB, etc.
- La Bella Beatle Bass sets can ONLY be used on Beatle basses (and similar models). Do NOT try them on standard basses (or else they will break).

Please email info@labella.com with specific questions about your bass!

Thanks for playing La Bella!

